

About installation of SGEMCS

- How to install

There are an object for 32 bit and an object for 64 bit in an installation file. Please use it according to the environment of the personal computer to install.

The installation file for 32 bit is contained into a "32-bit" folder, and the installation file for 64 bit is contained in the "64-bit" folder. Please perform "setup.exe" in each folder.

In the case of Windows Vista and 7, please install by "administrator authority."

If a user do not have an administrator authority, right-click setup.exe with a mouse. Then a pop up menu is displayed, click "performing as an administrator" and install.

- Version up

When you perform upgrade from the version earlier than Ver.6.0.0, please install in the following procedure.

1. Uninstall an older version.
2. Delete **ProInfo.dll** in the folder which the older version was installed
(usually C:\Program Files\SIGMA-KOKI\SGEMCS)
3. Install Ver.6.0.0. or later.

- License

When upgrade is performed from the version earlier than Ver.6.0.0, you need to perform license registration again.

Please demand the key (password) to open the lock from our company (an E-mail, FAX, Web).

- Execution

In order to perform the above-mentioned software by Windows Vista or 7, please perform by "administrator authority."

※If you do not perform by "administrator authority", neither data nor a movement environment is saved.