

\*\*\*\*\*

Software for automatic positioning and measurement  
About SGADVANC Ver2.13

Jun 10. 2021 SIGMA KOKI CO.,LTD. Technology Center

\*\*\*\*\*

## Environment

### ◆Hardware

Personal computers with Windows 8/8.1/10 (64bit)

CPU	: Core i3 or more
Memory	: 4GB or more (Recommend 8GB or more)
HardDisk	: More than 200MB free space
Drive	: CD-ROM(Installation possible drive)

### ◆Software

OS	: Microsoft Windows 8/8.1/10 (64bit)
Application	: Microsoft Excel2010 or later

### ◆Interface

Please pay attention because there is restriction by the environment and devices etc. where it is used, although it is supporting the interfacing of RS232C/GPIB/USB/LAN in this product.

RS232C	: Standard correspondence.
GPIB	: GPIB is only supported for an interface board from National Instruments. Before using the GPIB interface, the device driver for the GPIB card must be installed.
USB	: USB control is supported for SHOT series only. Before using the USB interface, the device

driver for the SHOT series must be installed.  
USB can be used in the personal computers with  
Windows 8/8.1/10 (64bit).

LAN : Standard correspondence.

◆Others

We are doing action confirmation in the following environment.

Excel2007/2010/2013.

We are doing action confirmation by the following stage controller  
about the USB interface.

SHOT-302GS/304GS, FINE-01gamma, HIT-M series

The guarantee and also inquiry regarding the action of this product  
ask to cause in above environment limited.